

RAPTOR Tests – quick info

SD Biosensor	BD Veritor	PHE virology
Current COVID symptoms	Within 5 days of symptom onset	
Nasopharyngeal swab	Nasal swab	Throat swab
Read results on test device between 15-30 mins after sample has been added	Read results on BD Veritor machine 15 mins after sample has been added	Swabs sent to PHE Colindale lab
Sample must be tested within 1 hour (if kept at room temp)	Swab must be processed (extracted into buffer) within 1 hour Once processed, sample must be added to test device within 30 mins	
Qualitative test	Qualitative test (machine read)	