

## RAPTOR Tests – quick info

Lumira Dx	UKHSA virology
<p><b>Current</b> COVID-like or flu-like symptoms, within 12 days of symptom onset</p>	
<p>Nasal swab</p>	<p>Nose and throat swab</p>
<p>Read results on test device between <b>10-15 mins</b> after sample has been added</p>	<p>Swabs sent to UKHSA Colindale lab in London</p>
<p>Sample can be kept for up to <b>5 hours</b> before being tested</p> <p>Please use the Lumira test strip within <b>15 mins</b> of opening</p> <p>When running the test on the Lumira device, please follow the screen instructions</p>	<p>Please post on the same day of testing</p>
<p>Qualitative test</p>	